

Anti-Phospho-CHUK-T23 antibody (STJ11101104)

STJ11101104

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1147
Gene Symbol	CHUK
Uniprot ID	IKKA_HUMAN
Immunogen	LGTGG
Immunogen Region	
Specificity	A synthetic phosphorylated peptide around T23 of human IKK Alpha (NP_001269.3).
Immunogen Sequence	LGTGG