

Anti-STAT3 antibody (130-715) (STJ11101195)

STJ11101195

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6774
Gene Symbol	STAT3
Uniprot ID	STAT3_HUMAN
Immunogen	
Immunogen Region	130-715
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 130-770 of human STAT3 (NP_644805.1). NHPTAAVTEKQQMLEQLQ DVRKRVQDLEQMKVVENLQ DDFDFNYKTLKSQGDMQDLN GNNQSVTRQKMQQLEQMLTA LDQMRRSIVSELAGLLSAME YVQKTLTDEELADWKRQIQ ACIGGPPNICLDRLNWLTS LAESQLQTRQIKKLEELQQ KVSYKGDPIVQHRPMLERI VELFRNLKMSAFVVERQPCM PMHPDRPLVIKTVGVQFTTKV RLLVKFPELNYQLKIKVCI
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081