

Anti-TFE3 antibody (191-440) [S0MR] (STJ11101240)
STJ11101240

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Clone ID	S0MR
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:4000
Range	IHC-P:1:100-1:400
	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	7030
Gene Symbol	TFE3
Uniprot ID	TFE3_HUMAN
Immunogen	
Immunogen Region	191-440
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 191-440 of human TFE3 (P19532).
Immunogen Sequence	LQQARRQQVKQYLSTTLGPK LASQALTPPPGPASQAQLPA PEAHHTTGPTGSAPNSPMAL LTIGSSSEKEIDDVIDEIS LESSYNDEMLSYLPGGTTGL QLPSTLPVSGNLLDVYSSQG VATPAITVSNSCPALPNIK REISETAKALLKERQKQDN HNLIERRRRFNINDRIKELG TLIPKSSDPEMRWNKGITLK ASVDYIRKLQKEQQRSKDLE SRQRSLEQANRSLQLRIQE



Immunohistochemistry analysis of TFE3 in paraffin-embedded Human acinar soft tissue sarcoma using TFE3 Rabbit monoclonal antibody (STJ11101240) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM Tris/EDTA buffer pH 9.0 before commencing with immunohistochemistry staining protocol.



Western blot analysis of lysates from HeLa cells, using TFE3 Rabbit monoclonal antibody (STJ11101240) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Western blot analysis of various lysates using TFE3 Rabbit monoclonal antibody (STJ11101240) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081