

Anti-CHMP3 antibody (65-220) [S3MR] (STJ11101263)

STJ11101263

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S3MR
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	100526767/51652
Gene Symbol	CHMP3
Uniprot ID	CHMP3_HUMAN
Immunogen	
Immunogen Region	65-220
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 65-220 of human VPS24 (Q9Y3E7).
Immunogen Sequence	KEMIRSRKAVSKLYASKAHM NSVLMGMKNQLAVLRVAGSL QKSTVMKAMQSLVKIPEIQ ATMRELSKEMMKAGIIIEML EDTFESMDDQEEMEEEAEME IDRILFEITAGALGKAPSKV TDALPEPEPPGAMAASEDEE EEEEEALEAMQSRLATL



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081