

## Anti-ULK1/ULK2 antibody (1-283) (STJ11101479)

STJ11101479

### GENERAL INFORMATION


|                          |                    |
|--------------------------|--------------------|
| <b>Product Type</b>      | Primary antibodies |
| <b>Short Description</b> |                    |
| <b>Applications</b>      | WB/ELISA           |
| <b>Host/Source</b>       | Rabbit             |
| <b>Reactivity</b>        | Human              |

### PRODUCT PROPERTIES

|                      |                                                                                                                                     |
|----------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>     | Polyclonal                                                                                                                          |
| <b>Clone ID</b>      |                                                                                                                                     |
| <b>Concentration</b> | Lot specific                                                                                                                        |
| <b>Conjugation</b>   | Unconjugated                                                                                                                        |
| <b>Purification</b>  | Affinity purification                                                                                                               |
| <b>Dilution</b>      | WB:1:500-1:2000                                                                                                                     |
| <b>Range</b>         | ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements. |
| <b>Formulation</b>   | PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.                                                                                    |
| <b>Isotype</b>       | IgG                                                                                                                                 |
| <b>Storage</b>       | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                      |
| <b>Instruction</b>   |                                                                                                                                     |

### TARGET INFORMATION

|                           |                                                                                                                                                                                                                                                                |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | <a href="#">8408</a><br><a href="#">9706</a>                                                                                                                                                                                                                   |
| <b>Gene Symbol</b>        | <a href="#">ULK1</a><br><a href="#">ULK2</a>                                                                                                                                                                                                                   |
| <b>Uniprot ID</b>         | <a href="#">ULK1_HUMAN</a><br><a href="#">ULK2_HUMAN</a>                                                                                                                                                                                                       |
| <b>Immunogen</b>          |                                                                                                                                                                                                                                                                |
| <b>Immunogen Region</b>   | 1-283                                                                                                                                                                                                                                                          |
| <b>Specificity</b>        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-283 of human ULK1/ULK2 (NP_003556.1).                                                                                                                                          |
| <b>Immunogen Sequence</b> | MEPGRGGTETVGKFEFSRKD LIGHGAFVVFVKGRHREKHD LEVAVKCINKKNLAKSQTLL GKEIKILKELKHENIVALYD<br>FQEMANSVYLVMEYCNGGDL ADYLHAMRTLSEDTIRLFLQ QIAGAMRLHSHKGIHRDLK PQNILLSNPAGRANPNSIR<br>VKIADFGFARYLQSNMMAAT LCGSPMYMAPEVIMSQHYDG KADLWSIGTIVYQCLTGKAP FQASSPQDLRLFYEKNKTL |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081