

Anti-FARP1 antibody (350-550) (STJ11101524)

STJ11101524

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	10160
Gene Symbol	FARP1
Uniprot ID	FARP1_HUMAN
Immunogen	
Immunogen Region	350-550
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 350-550 of human FARP1 (NP_005757.1).
Immunogen Sequence	YVKEGGHKVQFERKHSKIH SIRSLASQPTELNSEVLEQS QQSTSLTFGEGAESPGGQSC RRGKEPKVSAGEPGSHPSA PRRSPAGNKQADGAASAPTE EEEEWKDRQQSKPQQPPQ STGSLTGSPHSELVNSQG GVAPANVTLPNLSPDTKQA SPLISPLLNDQACPRTDDED EGRRKRFPD KAYFIAKEVS T