

## Anti-FAM160B2 antibody (74-303) (STJ11101525)

STJ11101525

### GENERAL INFORMATION





|                          |                    |
|--------------------------|--------------------|
| <b>Product Type</b>      | Primary antibodies |
| <b>Short Description</b> |                    |
| <b>Applications</b>      | WB/IF/ICC/IP/ELISA |
| <b>Host/Source</b>       | Rabbit             |
| <b>Reactivity</b>        | Human/Mouse/Rat    |

### PRODUCT PROPERTIES

|                      |                                                                                                                                    |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>     | Polyclonal                                                                                                                         |
| <b>Clone ID</b>      |                                                                                                                                    |
| <b>Concentration</b> | Lot specific                                                                                                                       |
| <b>Conjugation</b>   | Unconjugated                                                                                                                       |
| <b>Purification</b>  | Affinity purification                                                                                                              |
| <b>Dilution</b>      | WB:1:500-1:2000                                                                                                                    |
| <b>Range</b>         | IF/CC:1:50-1:200                                                                                                                   |
|                      | IP:0.5 Mu g-4 Mu g antibody for 400 Mu g-600 Mu g extracts of whole cells                                                          |
|                      | ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements |
| <b>Formulation</b>   | PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.                                                                                   |
| <b>Isotype</b>       | IgG                                                                                                                                |
| <b>Storage</b>       | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                     |
| <b>Instruction</b>   |                                                                                                                                    |

### TARGET INFORMATION

|                           |                                                                                                                                                                                                                                                 |
|---------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | 64760                                                                                                                                                                                                                                           |
| <b>Gene Symbol</b>        | FHIP2B                                                                                                                                                                                                                                          |
| <b>Uniprot ID</b>         | FHI2B_HUMAN                                                                                                                                                                                                                                     |
| <b>Immunogen</b>          |                                                                                                                                                                                                                                                 |
| <b>Immunogen Region</b>   | 74-303                                                                                                                                                                                                                                          |
| <b>Specificity</b>        | Recombinant fusion protein containing a sequence corresponding to amino acids 74-303 of human FAM160B2 (NP_073586.5).                                                                                                                           |
| <b>Immunogen Sequence</b> | GEAGPCLEYLLQHKILETLC TLGKAEYPPGMRQQVFQFFS KVLAQVQHPLLHYLSVHRPV QKLLRLGGTASGVSVEKEEV QFTTVLCSKIQQDPELLAYI LEGKKIVGRKKACGEPTALP KDTTSHGDKDCSHDGAPARP QLDGESCQAQALNSHMPAET EELDGGTTESNLITSLGLC QSKKSRVALKAQENLLLLVS MASPAATYLVQSSACCPAI VRHLCQLYRS |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081