

Anti-CCNA2 antibody (11-210) [S4MR] (STJ11101574)
STJ11101574

GENERAL INFORMATION





Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/IP/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S4MR
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:6000
Range	IHC-P:1:200-1:2000 IF/ICC:1:100-1:1000 IP:0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your s PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Formulation	
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	890
Gene Symbol	CCNA2
Uniprot ID	CCNA2_HUMAN
Immunogen	
Immunogen Region	11-210
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 11-210 of human Cyclin A2 (P20248).
Immunogen Sequence	TREAGSALLALQQTALQEDQ ENINPEKAAPVQQPRTRAAL AVLKSGNPRGLAQQRPKTR RVAPLKDLPVNDEHVTVPW KANSKQPAFTIHVDEAEKEA QKKPAESQKIEREDALAFNS AISLPGPRKPLVPLDYPMDG SFESPHTMDMSIILEDEKPV SVNEVPDYHEDIHTYLREME VKCKPKVGYMKKQPDITNSM



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081