

## Anti-NFKB2 antibody (745-877) [S5MR] (STJ11101725)

STJ11101725

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	
<b>Applications</b>	WB/IHC-P/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse/Rat

### PRODUCT PROPERTIES


<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	S5MR
<b>Concentration</b>	Lot specific
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution</b>	WB:1:1000-1:6000
<b>Range</b>	IHC-P:1:100-1:500
	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
<b>Formulation</b>	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	4791
<b>Gene Symbol</b>	NFKB2
<b>Uniprot ID</b>	NFKB2_HUMAN
<b>Immunogen</b>	
<b>Immunogen Region</b>	745-877
<b>Specificity</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 745-877 of human NF-Kappa B2 (Q00653).
<b>Immunogen Sequence</b>	LLLNAAQNTMEPPLTPPSPA GPGLSLGDTALQNLEQLLDG PEAQGSWAELAERLGLRSLV DTYRQTTSPSGSLLRSYELA GGDLAGLLEALSDMGLEEGV RLLRGPETRDKLPSTAIEVKE DSAYGSQSVEQEA



Immunohistochemistry analysis of NF-Kappa B2 in paraffin-embedded human appendix using NF-Kappa B2 Rabbit monoclonal antibody (STJ11101725) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Western blot analysis of various lysates using NF-Kappa B2 Rabbit monoclonal antibody (STJ11101725) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081