

Anti-BRD4 antibody (600-722) (STJ11101789)

STJ11101789

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23476
Gene Symbol	BRD4
Uniprot ID	BRD4_HUMAN
Immunogen	
Immunogen Region	600-722
Specificity	A synthetic peptide corresponding to a sequence within amino acids 600-722 of human BRD4 (NP_055114.1). ESEEDKCKPMSYEEKRQLS LDINKLPGEKLGRVWHIIQS REPSLKNSNPDEIEIDFETL KPSTLRELERYVTSCLRKKR
Immunogen Sequence	KPQAEKVDVIAGSSKMKGFS SSESESSSESSSSDSEDSET GPA



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081