

Anti-SLIT2 antibody (1400-1500) [S7MR] (STJ11101877)
STJ11101877

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Clone ID	S7MR
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:2000 IHC-P:1:200-1:800 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9353
Gene Symbol	SLIT2
Uniprot ID	SLIT2_HUMAN
Immunogen	
Immunogen Region	1400-1500
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1400-1500 of human Slit2 (O94813).
Immunogen Sequence	GHGGVLCDEEDLFNPCQAI KCKHKGKRLSGLGQPYCECS SGYTGDSCDREISCRGERIR DYYQKQQGYAACQTTKKVSR LECRGGCAGGQCCGPLRSKR R



Immunohistochemistry analysis of paraffin-embedded rat ovary using Slit2 Rabbit monoclonal antibody (STJ11101877) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM Tris/EDTA buffer pH 9.0 before commencing with immunohistochemistry staining protocol.



Western blot analysis of extracts of various cell lines, using Slit2 Rabbit monoclonal antibody (STJ11101877) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 3min.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081