

Anti-MYH11 antibody (1160-1313) [S8MR] (STJ11101988)
STJ11101988

GENERAL INFORMATION





Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S8MR
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IHC-P:1:500-1:5000 IF/ICC:1:400-1:1600 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	4629
Gene Symbol	MYH11
Uniprot ID	MYH11_HUMAN
Immunogen	
Immunogen Region	1160-1313
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1160-1313 of human SMMHC/MYH11 (NP_002465.1).
Immunogen Sequence	DSTATQQLRAKREQEVTVL KKALDEETRSHEAQVQEMRQ KHAQAVEELTEQLEQFKRAK ANLDKNKQTLEKENADLAGE LRVLGQAKQVEHKKKKLEA QVQELQSKCSDGERARAELN DKVHKLQNEVESVTGMLNEA EGKAIKLAKDVASL



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081