

Anti-BCL10 antibody (1-233) [S3MR] (STJ11102093)

STJ11102093

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S3MR
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	8915
Gene Symbol	BCL10
Uniprot ID	BCL10_HUMAN
Immunogen	
Immunogen Region	1-233
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-233 of human Bcl10 (O95999). MEPTAPSLTEEDLTEVKKDA LENLRVYLCEKIIAERHFDH LRAKKILSREDTEEISCRYS SRKRAGKLLDYLQENPKGLD
Immunogen Sequence	TLVESIRREKTNFLIQKIT DEVLKLRNIKLEHLKGLKCS SCEPPDGATNNLSRSNSDE SNFSEKLRASTVMYHPAGES STTPFFSTNSSNLPLVLEVG RTENTIFSSTTLPRPGDPGA PPLPPDLQLEEEGTCANSSE MFLPLRSRTVSRQ



Western blot analysis of extracts of various cell lines, using Bcl10 Rabbit monoclonal antibody (STJ11102093) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081