

Anti-TNNT2 antibody (199-298) [S5MR] (STJ11102195)

STJ11102195

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S5MR
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC-P:1:50-1:200 IF/ICC:1:50-1:200
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7139
Gene Symbol	TNNT2
Uniprot ID	TNNT2_HUMAN
Immunogen	
Immunogen Region	199-298
Specificity	A synthetic peptide corresponding to a sequence within amino acids 199-298 of human Cardiac troponin T (TNNT2) (P45379).
Immunogen Sequence	QKQAQTERKSGKRQTEREKK KKILAERRKVLAIHLDNQ LREKAKELWQSIYNLEAEKF DLQEKFKQKYEINVLNRRI NDNQKVSCTRGRKAKVTGRWK



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081