

Anti-T antibody (266-435) [S3MR] (STJ11102243)
STJ11102243

GENERAL INFORMATION




Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S3MR
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 IF/CC:1:50-1:200
	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6862
Gene Symbol	TBXT
Uniprot ID	TBXT_HUMAN
Immunogen	
Immunogen Region	266-435
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 266-435 of human Brachyury (O15178).
Immunogen Sequence	GALSLPSTHSCDRYPTLRSH RSSPYSPYAHRRNNSPTYSD NSPACLSMLQSHDNWSSLGM PAHPMSMLPVSHNASPPTSSS QYPSLWVSNGAVTPGSQAA AVSNGLG AQFFRGSPAHYTP LTHPVSAPSSSGSPLYEGAA AATDIVDSQYDAAAQGRLIA SWTPVSPSPM



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081