

Anti-RUNX1 antibody (354-453) [S0MR] (STJ11102260)

STJ11102260

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S0MR
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	860 861 864 RUNX2 RUNX1 RUNX2_HUMAN RUNX1_HUMAN RUNX3_HUMAN
Immunogen	354-453
Immunogen Region	
Specificity	A synthetic peptide corresponding to a sequence within amino acids 354-453 of human RUNX1/2/3 (NP_001001890.1).
Immunogen Sequence	LPPPYPGSSQAQGGPFQASS PSYHLYYGASAGSYQFSMVG GERSPPRILPPCTNASTGSA LLNPSLPNQSDVVEAEGSHS NSPTNMAPSARLEEAVWRPY



Western blot analysis of extracts of various cell lines, using RUNX1/2/3 Rabbit monoclonal antibody (STJ11102260) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081