

**Anti-HP antibody (307-406) [S3MR] (STJ11102373)**  
STJ11102373

**GENERAL INFORMATION**


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	
<b>Applications</b>	WB/IHC-P/IF/ICC/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse/Rat

**PRODUCT PROPERTIES**


<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	S3MR
<b>Concentration</b>	Lot specific
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB:1:500-1:2000 IHC-P:1:50-1:200 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 $\mu$ g/mL. Please optimize the concentration based on your specific assay requirements.
<b>Formulation</b>	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**


<b>Gene ID</b>	3240
<b>Gene Symbol</b>	HP
<b>Uniprot ID</b>	HPT_HUMAN
<b>Immunogen</b>	
<b>Immunogen Region</b>	307-406
<b>Specificity</b>	A synthetic peptide corresponding to a sequence within amino acids 307-406 of human Haptoglobin (HP) (P00738).
<b>Immunogen Sequence</b>	DQCI RHYEGSTVPEK KTPKS PVGVQPILNEHTFCAGMSKY QEDTCYGDAGSAFAVHDL EE DTWYATGILSFDKSCA VA EYGVYVKVTSIQDWVQKTI AEN



Western blot analysis of extracts of various cell lines, using Haptoglobin (HP) (HP) rabbit monoclonal antibody (STJ11102373) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 10s.



Immunofluorescence analysis of rat liver using Haptoglobin (HP) (HP) Rabbit monoclonal antibody (STJ11102373) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Immunohistochemistry of paraffin-embedded human liver cancer using Haptoglobin (HP) (HP) rabbit monoclonal antibody (STJ11102373) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM Tris/EDTA buffer pH 9.0 before commencing with immunohistochemistry staining protocol.

Immunohistochemistry analysis of paraffin-embedded human liver cancer using Haptoglobin (HP) (HP) Rabbit monoclonal antibody (STJ11102373) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM Tris/EDTA buffer pH 9.0 before commencing with immunohistochemistry staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081