

Anti-SOCS2 antibody (1-198) (STJ11102834)

STJ11102834

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	8835
Gene Symbol	SOCS2
Uniprot ID	SOCS2_HUMAN
Immunogen	
Immunogen Region	1-198
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-198 of human SOCS2 (NP_001257400.1).
Immunogen Sequence	MTLRCLPSGNGGEGTRSQW GTAGSAEESPQAARLAKAL RELGQTGWYWGSMTVNEAKE KLKEAPEGTFILRDSHSDY LLTISVKTSAGPTNLRIEQ DGKFRLDISIICVKSKLKQFD SVVHLIDYYVQMCKDKRTGP EAPRNGTVHLYLTKPLYTSA PSLQHLCLRTINKCTGAIWG LPLPTRLKDYLEEYKFQV