

Anti-Bcl9 antibody (1-500 aa) (STJ11102880)

STJ11102880

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IHC-P:1:50-1:200 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	77578
Gene Symbol	Bcl9
Uniprot ID	BCL9_MOUSE
Immunogen	
Immunogen Region	1-500 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-500 of mouse Bcl9. (Q9D219).
Immunogen Sequence	MHPSNPKVRSSPSGNTQSSP KSKQEVMPRPPTVMSPSGNP QLDSKFSNQGKPGGSASQSQ PSPCDSKSGGHTPKALPGPG GSMGLKNGAGNGAKGKGKRE RSISADSFQDRDPGTPNDDS DIKCNADHIKSQESQHTP HSMTPSTATAPRSSTPSHGQ TPAPEPISAQKTPAKVVYVF STEMANKAAEAVLKQGVETI VSFHIQNISNSKSERSTAPL NTQIPTLRNDPKPLPQQPP



Immunofluorescence analysis of U2OS cells using Bcl9 Rabbit polyclonal antibody (STJ11102880) at dilution of 100 (40x lens). Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081