

Anti-CC2D1A antibody (750-951) [S4MR] (STJ11103014)

STJ11103014

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S4MR
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IHC-P:1:50-1:200
	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54862
Gene Symbol	CC2D1A
Uniprot ID	C2D1A_HUMAN
Immunogen	
Immunogen Region	750-951
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 750-951 of human CC2D1A (Q6P1N0).
Immunogen Sequence	VLGTAQLKLDALAEACEVRE ILEVLDGRRPTGGRLEVMMVR IREPLTAQQLETTTERWLVI DPVPAAVPTQVAGPKGKAPP VPAPARESGNRSARPLHSLVLAQFDQERLERKILALRQAR RPVPPEVAQYQDIMQRSQW QRAQLEQGGVGIRREYAAQL ERQLQFYTEAARRLGNDGSR DAAKEALYRRNLVESELQRL RR



Immunohistochemistry analysis of paraffin-embedded human colon using CC2D1A Rabbit monoclonal antibody (STJ11103014) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM Tris/EDTA buffer pH 9.0 before commencing with immunohistochemistry staining protocol.