

Anti-NAT10 antibody (750-1025) [S6MR] (STJ11103016)
STJ11103016

GENERAL INFORMATION





Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S6MR
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:4000
Range	IHC-P:1:200-1:800 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55226
Gene Symbol	NAT10
Uniprot ID	NAT10_HUMAN
Immunogen	
Immunogen Region	750-1025
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 750-1025 of human NAT10 (Q9H0A0).
Immunogen Sequence	KLTDEDEADQGGWLAFAWK DFRRLFLALLSYQFSTFSPS LALNIIQNRNMGKPAQPALS REELEALFLPYDLKRLEMYS RNMVDYHLIMDMIPASRIY FLNQLGDLALSAAQSALLG IGLQHKSVQDQLEKEIELPSG QLMGLFNRIIRKVVKLFNEV QEKAIEEQMVAAKDVVMEPT MKTLSDDLDEAAKEFQEKHK KEVGKLKSMDDLSEYIIRGDD EEWNEVLNKAGPNASIISL



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081