

Anti-ANKRD27 antibody (560-710) (STJ11103248)

STJ11103248

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84079
Gene Symbol	ANKRD27
Uniprot ID	ANR27_HUMAN
Immunogen	
Immunogen Region	560-710
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 560-710 of human ANKRD27 (NP_115515.2).
Immunogen Sequence	IGNEKGDTPHIAARWGYQG VIETLLQNGASTEIQNRLKE TPLKCALNSKILSVMEAYHL SFERRQKSSEAPVQSPQRSV DSISQESSTSSFSSMSASSR QEETKKDYREVEKLLRAVAD GDLEMVRYLLEWTEEDLEDA EDTVSAADPEF