

Anti-AMER1 antibody (550-660) [S7MR] (STJ11103307)

STJ11103307

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S7MR
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	139285
Gene Symbol	AMER1
Uniprot ID	AMER1_HUMAN
Immunogen	
Immunogen Region	550-660
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 550-660 of human AMER1 (NP_689637.3).
Immunogen Sequence	PPGAMETEEERLVTIQKQLL YWELRREQLEAQEARAREAH AREAHAREAYTREAYGREAY AREAHTWEAHGREARTREAY AREVRCRETQVRETQARQEK PVLEYQMRPLG