

Anti-TET3 antibody (1-285) (STJ11103351)

STJ11103351

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	194388
Gene Symbol	Tet3
Uniprot ID	TET3_MOUSE
Immunogen	
Immunogen Region	1-285
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-285 of mouse TET3 (Q8BG87).
Immunogen Sequence	MSQFQVPLAVQPDLSGLYDF PQGQVMVGGFQGPGLPMAGS ETQLRGGGDGRKKRRCGTC DPCRRLENCGSCTSCNRRRT HQICKLRKCEVLKKKAGLLK EVEINAREGTGPWAQGATVK TGSELSPVDGPVPGQMDSGP VYHGDSRQLSTSGAPVNGAR EPAGPGLLGAAGPWRVDQKP DWEAASGPTHAAARLEDAHDL VAFSAVAEAVSSYGALSTRL YETFNREMSREAGSNRGP



Western blot analysis of lysates from wild type (WT) and 293T cells transfected with TET3 using TET3 Rabbit polyclonal antibody (STJ11103351) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ0000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081