

Anti-Desmoglein 2 antibody (1001-1118) [S6MR] (STJ11103566)

STJ11103566

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S6MR
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1829
Gene Symbol	DSG2
Uniprot ID	DSG2_HUMAN
Immunogen	
Immunogen Region	1001-1118
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1001-1118 of human Desmoglein 2 (Q14126).
Immunogen Sequence	PLEGTQHLQDVPYVMVRERE SFLAPSSGVQPTLAMPNIAV GQNVTVTERVLAPASTLQSS YQIPTENSM TARNTTVSGAG VPGPLPDFGLEESGHSNSTI TTSSTRVTKHSTVQHSYS



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081