

Anti-Beta-Tubulin antibody (40-260) (STJ11103688)

STJ11103688

GENERAL INFORMATION




Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB: 1:500-1:2000 IHC-P: 1:50-1:200 IF/ICC: 1:50-1:200
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	203068
Gene Symbol	TUBB
Uniprot ID	TBB5_HUMAN
Immunogen	
Immunogen Region	40-260
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 40-260 of human Beta-Tubulin (NP_006077.2).
Immunogen Sequence	SDLQLERISVYYNEASSHKY VPRAILVDLEPGTMDSVRSG AFGHLFRPDNFIQGSGAGN NWAKGHYTEGAELVDSVLDV VRKECENCDCQLQGFLTHSL GGGTGSGMGTLLISKVREEY PDRIMNTFSVWPSPKVSDTV VEPYNATLSIHQLVENTDET YCIDNEALYDICFRTLKLAT PTYGDLNHLVSATMSGVTTSLRFPGLNADLRKLAVNMVP F



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081