

## Anti-GST-Tag antibody (1-239) [S0MR] (STJ11103710)

STJ11103710

### GENERAL INFORMATION



<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	
<b>Applications</b>	WB/IP/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Species independent

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	S0MR
<b>Concentration</b>	Lot specific
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution</b>	WB:1:2000-1:6000
<b>Range</b>	IP:3 Mu g antibody for 1 Mu g extracts of recombinant protein ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
<b>Formulation</b>	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	
<b>Immunogen Region</b>	1-239
<b>Specificity</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 1-218 of schja GST protein (P08515).
<b>Immunogen Sequence</b>	MSPILGYWKIKGLVQPTRLL LEYLEEKYEELHYERDEGDK WRNKKFELGLEFPNLPYYID GDVKLTQSMAIRYIADKHN MLGGCPKERAIEISMLEGAVL DIRYGVSR IAYSKDFETLKV DFLSKLPEMLKMFEDRLCHK TYLNGDHVTHPDFMLYDALD VLYMDPMCLDAFPKLVCFK KRIEAIQIDKYLKSSKYIA WPLQGWQATFGGGDHPPK



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081