

Anti-Phospho-Tyk2-Y292 antibody (STJ11103717)

STJ11103717

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7297
Gene Symbol	TYK2
Uniprot ID	TYK2_HUMAN
Immunogen	PCYIR
Immunogen Region	
Specificity	A synthetic phosphorylated peptide around Y292 of human Tyk2 (NP_003322.3).
Immunogen Sequence	PCYIR