

Anti-TOP2A antibody (1314-1503) (STJ11103810)

STJ11103810

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7153
Gene Symbol	TOP2A
Uniprot ID	TOP2A_HUMAN
Immunogen	
Immunogen Region	1314-1503
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1314-1503 of human DNA topoisomerase II alpha (TOP2A) (NP_001058.2).
Immunogen Sequence	ETEPRAATKTKFTMDLSD EDFSDFDEKTDDEDFVPSDA SPPKTKTSPKLSNKLKPQK SVVSDLEADDVKGSVPLSSS PPATHFPDETEITNPVKKK VTKKTAAKSQSSTSTTGAK KRAAPKGTKRDPALNSGVQS KPDPAKTKNRRKRKPSTSD SDSNFEKIVSKAVTSKSKG ESDDFHMDFD



Immunofluorescence analysis of U-2 OS cells using DNA topoisomerase II alpha (DNA topoisomerase II alpha (TOP2A)) Rabbit polyclonal antibody (STJ11103810) at dilution of 1:100 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081