

Anti-OsELF3 antibody (1-172 aa) (STJ11103914)
STJ11103914

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	O.sativa

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen	1-172 aa
Region	
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-172 of Oryza sativa OsELF3. (Q9SNQ6).
Immunogen	MATRGGGGGGGKEAKGKVM GPLFPRLHVNDAAKGGGPRA PPRNKMALYEQFTVPSHRFS GGGGGGGVGGSPAHSSTAAS
Sequence	QSQSQSQVYGRDSSLFPFN VPSNRPGHSTEKINSKINK KISGSRKELGMLSSQTKGMD IYASRSTAEAPQRRRAENTIK SSSGKRLADDE



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081