

Anti-VCAN antibody (1400-1586) (STJ11103954)

STJ11103954

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1462
Gene Symbol	VCAN
Uniprot ID	CSPG2_HUMAN
Immunogen	
Immunogen Region	1400-1586
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1400-1586 of human Versican (NP_004376.2).
Immunogen Sequence	TDVTTTPSVQYINGKHLVTT VPKDPEAAEARRGQFESVAP SQNFSDSSESDTHPFVIAKT ELSTAVQPNESTETTESLEV TWKPETYPETSEHFSGGEPD VFPTVPFHEEFESGTAKKGA ESVTERDTEVGHQAHEHTEP VSLFPEESSGEIAIDQESQK IAFARATEVTFGEEVEKSTS VTYTPTI



Immunofluorescence analysis of paraffin-embedded mouse brain using Versican Rabbit polyclonal antibody (STJ11103954) at dilution of 1:100 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of paraffin-embedded rat brain using Versican Rabbit polyclonal antibody (STJ11103954) at dilution of 1:100 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081