

Anti-TNKS1BP1 antibody (700-1000) (STJ11103972)

STJ11103972

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IHC-P:1:50-1:200
	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	85456
Gene Symbol	TNKS1BP1
Uniprot ID	TB182_HUMAN
Immunogen	
Immunogen Region	700-1000
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 700-1000 of human TNKS1BP1 (NP_203754.2).
Immunogen Sequence	RRGAGAEKDTQSPSTCSEG LLGWSQKDLQSEFGITGDPQ PSSFSPPSWCQGASQDYGLG GASPRGDPGLGERDWTSKYG QGAGEGSTREWASRCGIGQE EMEASSSQDQSKVSAPGVLTAQDRVVGKPAQLGTQRSQEA DVQDWEFRKRDSQGTYSRRD AELQDQEFGKRDSLGTYSRR DVS LGDWEFGKRDSL GAYAS QDANEQQGDLGKRDDHHGRYS SQDADEQDWEFQKRDVSLG



Immunohistochemistry analysis of paraffin-embedded rat stomach using TNKS1BP1 Rabbit polyclonal antibody (STJ11103972) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded mouse lung using TNKS1BP1 Rabbit polyclonal antibody (STJ11103972) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded human colon using TNKS1BP1 Rabbit polyclonal antibody (STJ11103972) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081