

Anti-ZNF213 antibody (132-459) (STJ11104059)
STJ11104059

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	7760
Gene Symbol	ZNF213
Uniprot ID	ZN213_HUMAN
Immunogen	
Immunogen Region	132-459
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 132-459 of human ZNF213 (NP_004211.1).
Immunogen Sequence	RQDVPSSEEAPEAAGRGSA TGPPTVGARRRPSVPQEQH SHSAQPPALLKEGRPGETTD TCFVSGVHGPVALGDIPFYF SREEWGTLDPAQRDLFWDIK RENSNTTLGFGGLKGQSEKS LLQEMVPVVPGGTGSDVTVS WSPPEAEAWESENRPRAALG PVVGARRGRPPTRRRQFRDL AAEKPHSCGQCGRFRWGS D LARHQRTHTGKPHKCECD KSFRSSSDLVRHQGVHTGE




Immunohistochemistry analysis of paraffin-embedded rat testis using ZNF213 Rabbit polyclonal antibody (STJ11104059) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded mouse brain using ZNF213 Rabbit polyclonal antibody (STJ11104059) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded human lung cancer using ZNF213 Rabbit polyclonal antibody (STJ11104059) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma using ZNF213 Rabbit polyclonal antibody (STJ11104059) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081