

Anti-ZHX2 antibody (260-500) [S4091RM] (STJ11104091)

STJ11104091

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S4091RM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	22882
Gene Symbol	ZHX2
Uniprot ID	ZHX2_HUMAN
Immunogen	
Immunogen Region	260-500
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 260-500 of human ZHX2 (Q9Y6X8). VPLNNTTKYNSALDTNATMIN SFNKFYPYPTQAELSWLTAAS KHPEEHIRIWFATQRLKHGI SWSPEEVEEARKKMFNGTIQ SVPPITITVLPAPLAPTKVTQ PILQTALPCQILGQTSVLVT QVTSGSTTVSCSPITLAVAG VTNHGGKRPLVTPQAAPEPK RPHIAQVPEPPPKVANPPLT PASDRKKTKEQIAHLKASFL QSQFPDDAEVYRLIEVTGLA RSEIKKWFS DhRYRCQRGI
Sequence	



Western blot analysis of various lysates, using ZHX2 antibody (STJ11104091) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 45s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081