

Anti-CRTC2 antibody (594-693) [S4099RM] (STJ11104099)

STJ11104099

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S4099RM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	200186
Gene Symbol	CRTC2
Uniprot ID	CRTC2_HUMAN
Immunogen	
Immunogen	594-693
Region	
Specificity	A synthetic peptide corresponding to a sequence within amino acids 594-693 of human TORC2 (Q53ET0).
Immunogen	PQDPHTFNHQNLTHCSRHGS GPNILTDSSPGFSKEIAA ALAGVPGFEVSAAGLELGLG LEDELRLMEPLGLEGLNMLSD
Sequence	PCALLPDPAVEESFRSDRLQ



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081