

Anti-USF1 antibody (100-200) [S4102RM] (STJ11104102)

STJ11104102

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S4102RM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7391
Gene Symbol	USF1
Uniprot ID	USF1_HUMAN
Immunogen	
Region	100-200
Specificity	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human USF1 (P22415).
Immunogen Sequence	TSDDAVDTGTAETHYTYF PSTAVGDGAGGTTSGSTAAV VTTQGSALLGQATPPGTGQ FVMMSPQEVLGGSQRSIA PRTHPYSPKSEAPRTTRDEK R



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081