

Anti-Pseudouridine/5-ribosyluracil antibody [S4148RM] (STJ11104148)

STJ11104148

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	DB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Species independent

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S4148RM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	DB:1:500-1:1000 IHC-P:1:500-1:1000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	Chemical compounds corresponding to Pseudouridine/5-ribosyluracil/Y.
Immunogen Sequence	



Dot-blot analysis of different sorts of chemical compounds using Pseudouridine/ 5-ribosyluracil Rabbit monoclonal antibody (STJ11104148) at 1:1000 dilution.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081