

Anti-Ezrin/Radixin/Moesin antibody (487-586) (STJ11104178)
STJ11104178

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7430 5962 4478 EZR RDX
Uniprot ID	EZRI_HUMAN RADI_HUMAN MOES_HUMAN
Immunogen	
Immunogen Region	487-586
Specificity	A synthetic peptide corresponding to a sequence within amino acids 487-586 of human Ezrin/Radixin/Moesin (NP_003370.2). ESLQDEGAEPTGYSAELSSSE GIRDDRNEEKRITEAEKNER VQRQLLTLSSELSQARDENK RTHNDIIHNENMRQGRDKYK
Immunogen Sequence	TLRQIRQNTKQRIDEFEAL



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081