

Anti-PMF1 antibody (1-205) [S4199RM] (STJ11104199)
STJ11104199

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Clone ID	S4199RM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	100527963/11243
Gene Symbol	PMF1
Uniprot ID	PMF1_HUMAN
Immunogen	
Immunogen Region	1-205
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-205 of human PMF1 (NP_009152.2).
Immunogen Sequence	MAEASSANLGS GCEEKRHEG SSSESVPPTTISRVLKLD MVDFTLQKLVAAGSYQRFTD CYKCFYQLQPAMTQQIYDKF IAQLQTSIREEISDIKEEGN LEAVLNALDKIVEEGKVRKE PAWRPSGIPEKDLHSMAPY FLQQRDTLRRHVQKQEAENQ QLADAVLAGRRQVEELQLQV QAQQAWQALHREQRELVAV LREPE



Western blot analysis of various lysates using PMF1 Rabbit monoclonal antibody (STJ11104199) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 180s.



Western blot analysis of various lysates using PMF1 Rabbit monoclonal antibody (STJ11104199) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081