

Anti-COPB2 antibody (657-906) [S4264RM] (STJ11104264)

STJ11104264

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S4264RM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:3000-1:20000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9276
Gene Symbol	COPB2
Uniprot ID	COPB2_HUMAN
Immunogen	
Immunogen Region	657-906
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 657-906 of human COPB2 (NP_004757.1). AYQLAVEAESEKWKQLAEL AISKQFGLAQECLHHAQDY GGLLLLATASGNANMVNKL AEGAERDGGKNNVAFMSYFLQG KVDACLELLIRTGRLPEAAF LARTYLPQSVSRVVKLWREN LSKVNQAAESLADPTEYEN LFPGLKEAFVVEEWWKETHA DLWPAKQYPLVTPNEERNVM EEGKDFQPSRSTAQQLDGK PASPTPVIVASHTANKEEKS LLELEVLDNLELEDIDTT
Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081