

Anti-DSG2 antibody (635-935) (STJ11104316)

STJ11104316

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1829
Gene Symbol	DSG2
Uniprot ID	DSG2_HUMAN
Immunogen	
Immunogen Region	635-935
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 635-935 of human DSG2 (NP_001934.2).
Immunogen Sequence	CHCGKGAAGFTPIPGTIEML HPWNNEGAPPEDKVPSFLP VDQGGSLVGRNGVGGMAKEA TMKGSSSASIVKGQHEMSEM DGRWEEHRSLLSGRATQFTG ATGAIMTTETTKTARATGAS RDMAGAAQAVALNEEFLRN YFTDKAASYTEEDENHTAKD CLLVYSQEETESLNASIGCC SFIEGELDDRFLDDLGLKFK TLAEVCLGQKIDINKEIEQR QKPATETSMNTASHSLCEQ



Western blot analysis of various lysates using DSG2 Rabbit polyclonal antibody (STJ11104316) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081