

Anti-FLOT2 antibody (239-428) (STJ11104327)

STJ11104327

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2319
Gene Symbol	FLOT2
Uniprot ID	FLOT2_HUMAN
Immunogen	
Immunogen Region	239-428
Specificity	Recombinant Protein corresponding to a sequence within amino acids 239-428 of human FLOT2 (NP_004466.2). LAYELQGAREQQKIRQEEIE IEVQRKKQIAVEAQEILRT DKELIATVRRPAEAEAHRIQ QIAEGEKVKQVLLAQAEAEK
Immunogen Sequence	IRKIGEEAAVIEAMGKAEA ERMKLKAEAYQKYGDAAKMA LVLEALPQIAAKIAAPLTKV DEIWLSGDNSKVTSEVNRL LAELPASVHALTGVDLSKIP LIKKATGVQV



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081