

Anti-NR3C1 antibody (1-280) (STJ11104346)

STJ11104346

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	2908
Gene Symbol	NR3C1
Uniprot ID	GCR_HUMAN
Immunogen	
Immunogen Region	1-280
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human Glucocorticoid Receptor (NP_001191194.1).
Immunogen Sequence	MDSKESLTPGREENPSSVLA QERGDVDMFYKTLRGGATVK VSASSPSLAVASQSDSKQRR LLVDFPKGSVSNAQQPDLSK AVLSMGLYMGETETKVMGN DLGFPQQGQISLSSGETDLK LLEESIANLNRSTSVPENPK SSASTAVSAAPTEKEFPKTH SDVSSSEQHLKGQTGTNGGN VKLYTTDQSTFDILQDLEFS SGSPGKETNESPWRSDLLID ENCLLSPLAGEDDSFLLEG



Immunohistochemistry analysis of Glucocorticoid Receptor in paraffin-embedded rat liver using Glucocorticoid Receptor Rabbit polyclonal antibody (STJ11104346) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of Glucocorticoid Receptor in paraffin-embedded human tonsil using Glucocorticoid Receptor Rabbit polyclonal antibody (STJ11104346) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Western blot analysis of various lysates using Glucocorticoid Receptor Rabbit polyclonal antibody (STJ11104346) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081