

Anti-PBK antibody (1-322) (STJ11104418)

STJ11104418

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	55872
Gene Symbol	PBK
Uniprot ID	TOPK_HUMAN
Immunogen	
Immunogen Region	1-322
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-322 of human PBK (NP_060962.2).
Immunogen Sequence	MEGISNFKTPSKLSEKKKSV LCSTPTINIPASPFMQKLG GTGVNVYLMKRSPRGLSHSP WAVKKINPICNDHYRSVYQK RLMDEAKILSLHHPNIVGY RAFTEANDGSLCLAMEYGG KSLNDLIEERYKASQDPFPA AILKVALNMARGLKYLHQE KLLHGDIKSSNVVVGKDFE TIKICDVGVSPLDENMTVT DPEACYIGTEPWKPKEAVEE NGVITDKADIFAFGLTLWE



Western blot analysis of various lysates using PBK Rabbit polyclonal antibody (STJ11104418) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3 % nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1 s.



Western blot analysis of lysates from HeLa cells using PBK Rabbit polyclonal antibody (STJ11104418) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3 % nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1 s.



Western blot analysis of lysates from wild type (WT) and PBK knockout (KO) 293T cells using PBK Rabbit polyclonal antibody (STJ11104418) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3 % nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1 s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081