

Anti-CSNK1E antibody (247-416) (STJ11104560)

STJ11104560

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IHC-P: 1:50-1:200
Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	102800317/1454
Gene Symbol	CSNK1E
Uniprot ID	KC1E_HUMAN
Immunogen	
Immunogen Region	247-416
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 247-416 of human CSNK1E (NP_689407.1).
Immunogen Sequence	EFSTYLNFCRSLRFDDKPDY SYLRQLFRNLFHRQGFSDY VFDWNMLKFGAARNPEDVDR ERREHEREERMGLRGSATR ALPPGPPTGATANRLRSAAE PVASTPASRIQPAGNTSPRA ISRVDREKRVSMRLHRGAPA NVSSDLTGQRQEVSRIPASQ TSVPFDLHLGK



Immunohistochemistry analysis of paraffin-embedded rat ovary using CSNK1E Rabbit polyclonal antibody (STJ11104560) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded mouse kidney using CSNK1E Rabbit polyclonal antibody (STJ11104560) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded human liver using CSNK1E Rabbit polyclonal antibody (STJ11104560) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.