

Anti-PLIN2 antibody (300-437) (STJ11104660)

STJ11104660

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IF/CC:1:50-1:200
	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	123
Gene Symbol	PLIN2
Uniprot ID	PLIN2_HUMAN
Immunogen	
Immunogen Region	300-437
Specificity	A synthetic peptide corresponding to a sequence within amino acids 300-437 of human Perilipin-2 (NP_001113.2).
Immunogen Sequence	SHCAEHIESRTLAIARNLTQ QLQTTCHTLLSNIQGVQNI QDQAKHMGVMAGDIYSVFRN AASFKEVSDSLTSSKGQLQ KMKESLDDVMDYLVNNTPLN WLWGPFPQLTESQNAQDQG AEMDKSSQETQRSEHKTH



Immunofluorescence analysis of PC-12 cells using Perilipin-2 Rabbit polyclonal antibody (STJ11104660) at dilution of 1:50 (40x lens). Blue: DAPI for nuclear staining.



Western blot analysis of extracts of HuH-7, using Perilipin-2 Rabbit polyclonal antibody (STJ11104660) at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081