

Anti-IL10 antibody (19-178) [S4793RM] (STJ11104793)

STJ11104793

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S4793RM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:5000
Range	
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	25325
Gene Symbol	IL10
Uniprot ID	IL10_RAT
Immunogen	
Immunogen Region	19-178
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 19-178 of rat IL10 (NP_036986.2).
Immunogen Sequence	SKGHSIRGDNNCTHFPVSQT HMLRELRAAFSQVKTFFQKK DQLDNILLTDSLQDFKGYL GCQALSEMIKFYLVEVMPQA ENHGPEIKEHLNSLGEKLT LWIQLRRCHRFLPCENKSKA VEQVKNDFNKLQDKGVYKAM NEFDIFINCIEAYVTLKMKN



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081