

Anti-TOMM70 antibody (245-608) [S4820RM] (STJ11104820)

STJ11104820

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S4820RM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9868
Gene Symbol	TOMM70
Uniprot ID	TOM70_HUMAN
Immunogen	
Immunogen Region	245-608
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 245-608 of human TOM70 (NP_055635.3).
Immunogen Sequence	KNREPLMPSPQFIKSYFSSF TDDIISQPMLKGEKSDKDK KEGEALEVKENSGYLKAKQY MEEENYDKIIECSKEIDAE GKYMAEALLLRATFYLLIGN ANAAKPDLDKVISLKEANVK LRANALIKRGSMYMQQQQPL LSTQDFNMAADIDPQNADV HHRGQLKILLDQVEEAVADF DECIRLRPESALAAQAKCFA LYRQAYTGNNSSQIQAAMKG FEEVIKKFPRCAEGYALYA



Western blot analysis of various lysates, using TOMM70 antibody (STJ11104820) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081