

Anti-IPO5 antibody (20-300 aa) [S4950RM] (STJ11104950)

STJ11104950

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S4950RM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:2000-1:10000
Range	IF/CC:1:50-1:200
	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3843
Gene Symbol	IPO5
Uniprot ID	IPO5_HUMAN
Immunogen	
Immunogen Region	20-300 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 20-300 of human Karyopherin Beta 3 (IPO5). (NP_002262.4).
Immunogen Sequence	SPDNVVRKQAEETYENIPGQ SKITFLQAIRNTTAAEEAR QMAAVLLRRLSSAFDEVYP ALPSDVQTAIKSELLMIQM ETQSSMRKKVCDIAAELARN LIDEDGNNQWPEGLKFLFDS VSSQNVGLREAALHIFWNFP GIFGNQQQHYLDVIKRLVQ CMQDQEHPSIRTLSARATAA FILANEHNVALFKHFADLLP GFLQAVNDSCYQNDSDSVLKS LVEIADTVPKYLRPHLEAT



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081