

Anti-CD201 antibody (18-210 aa) (STJ11104990)

STJ11104990

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	19124
Gene Symbol	Procr
Uniprot ID	EPCR_MOUSE
Immunogen	
Immunogen Region	18-210 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 18-210 of mouse CD201. (NP_035301.2). LCNSDGSQSLHMLQISYFQD NHHVVRHQGNASLGKLLHTL EGPSQNVTLQLQPWQDPES WERTESGLQIYLTQFESLVK
Immunogen Sequence	LVYRERKENVFFPLTVSCSL GCELPEEEEEGSEPHVFFDV AVNGSAFVSFRPKTAVVWWSG SQEPSKAANFTLKQLNAYNR TRYELQEFLQDTCVEFLENH ITTQNMKGSTGR



Western blot analysis of Mouse spleen, using CD201 antibody (STJ11104990) at 1:600 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 180s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081